

## HEPATEK

### Halbautomatische Biopsie-Nadel, Menghini-Typ *Menghini type semi-automatic biopsy needle*

Automatische Aspiration. Mandrin mit Trokarspitze und am Spritzen-Kolben befestigt. Schneide-Spitze, Menghini-Typ, Inklusive Tiefenmarkierungen und Stopper. Innerer Ultraschall-Marker

*Automatic suction. Trocar tip stylet connected at the syringe plunger. Menghini type cutting needle tip. Depth marks and stopper included. Echogenic marker: inner*

VAKUUM-UNTERSTÜTZT



Trokar-Spitze  
*Trocar tip*



Innerer Ultraschall-Marker  
*Inner echogenic marker*

	Ø Gauge	Ø mm	Länge Length 40 mm	Länge Length 70 mm	Länge Length 90 mm	Länge Length 100 mm	Länge Length 120 mm	Länge Length 150 mm	Länge Length 200 mm
○	16	1,60	HS16-04	HS16-07	HS16-09	HS16-10	HS16-12	HS16-15	HS16-20
●	17	1,47	HS17-04	HS17-07	HS17-09	HS17-10	HS17-12	HS17-15	HS17-20
●	18	1,27	HS18-04	HS18-07	HS18-09	HS18-10	HS18-12	HS18-15	HS18-20
●	19	1,10	HS19-04	HS19-07	HS19-09	HS19-10	HS19-12	HS19-15	HS19-20
●	20	0,90	HS20-04	HS20-07	HS20-09	HS20-10	HS20-12	HS20-15	HS20-20
●	21	0,82	HS21-04	HS21-07	HS21-09	HS21-10	HS21-12	HS21-15	HS21-20
●	22	0,70	HS22-04	HS22-07	HS22-09	HS22-10	HS22-12	HS22-15	HS22-20

Box mit 10 Stück / Weitere Größen sind auf Anfrage lieferbar  
*Box 10 units / Different sizes are available on request*