

Available in 3 lengths

● Scalpel (36cm)

5mm laparoscopic Curved shear



● Scalpel (23cm)

5mm Curved shear for Open Procedures



● Scalpel (14cm)

5mm Curved shear for Open Procedures



Order Information and Specifications

Order Information

Product code	Product Description	Packaging
CSUS 6000	Generator with foot switch and test tip	1 pc/case
CS3605H/P	5mm Laparoscopic Curved Shear,Pistol, 36cm	3pc/case
CS2305H/P	5mm Open Curved Shear, Pistol, 23cm	3pc/case
CS1405H/P	5mm Open Curved Shear, Pistol, 14cm	3pc/case
TRA 6	Transducer	1 pc/case
TRASE/TRAS	Transducer cable	1 pc/case
CART 6000	Generator Cart	1 pc/case

Specifications

Generator dimensions	Generator Weight	Output Settings	frequency	Cart Dimensions
16cm(H)x36cm(w)x38cm(D)	10kg	5 Levels	55.5KHz	94cm(H)x39cm(W)x60cm(D)



Cart



Transducer



Foot Switch