

Alguide™

ENDOVASCULAR GUIDING CATHETER

The Alguide™ endovascular guiding catheter offers excellent torque and kink resistance, enabling fast and effective procedures.




- Outstanding torque control and kink resistance for an easier catheter positioning and handling
- Excellent curve retention and back-up support for challenging procedures
- Large inner lumen and wide variety of shapes to meet all requirements

TECHNICAL DATA

Diameters:	5F, 6F, 7F & 8F	Material main shaft:	Nylon
Lengths:	45, 60, 65, 80 & 100 cm	Braiding:	Stainless steel
Shaft:	5 segments (variable stiffness) with teflon inner liner	Tip material:	Pebax
		Radiopaque component:	Barium Sulfate

DIMENSIONS



ORDER INFORMATION ALGUIDE™

	Vertebral/ Multipurpose		Visceral/ Renal		Visceral/ Renal		Visceral/ Renal		Cerebral/ Visceral		Cerebral/ Visceral			
	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number		
5 Pack	5 F	100	503110202444	100	503110202462	100	503110202463	100	503110202464	100	503110202494	100	503110202495	
				80	503110202524	80	503110202526	80	503110202528	65	503110202515	65	503110202517	
				65	503110202523	65	503110202525	65	503110202527					
		6 F	100	503110202744	100	503110202762	100	503110202763	100	503110202764	100	503110202794	100	503110202795
					80	503110202824	80	503110202826	80	503110202828	65	503110202815	65	503110202817
					65	503110202823	65	503110202825	65	503110202827				
	7 F	100	503110203044	100	503110203062	100	503110203063	100	503110203064	100	503110203094	100	503110203095	
				80	503110203124	80	503110203126	80	503110203128	65	503110203115	65	503110203117	
				65	503110203123	65	503110203125	65	503110203127					
	8 F	100	503110203344	100	503110203362	100	503110203363	100	503110203364	100	503110203394	100	503110203395	
				80	503110203424	80	503110203426	80	503110203428	65	503110203415	65	503110203417	
				65	503110203423	65	503110203425	65	503110203427					

	Cerebral		Cerebral		Cerebral		Cerebral		Vertebral/ Multipurpose		Visceral/ Renal		
	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	
5 Pack	5 F	100	503110202484	100	503110202490	100	503110202485	100	503110202491	100	503110202474	100	503110202472
	6 F	100	503110202784	100	503110202790	100	503110202785	100	503110202791	100	503110202774	100	503110202772
	7 F	100	503110203084	100	503110203090	100	503110203085	100	503110203091	100	503110203074	100	503110203072
	8 F	100	503110203384	100	503110203390	100	503110203385	100	503110203391	100	503110203374	100	503110203372

ORDER INFORMATION ALGUIDE™ (continued)

	Visceral / Renal		Visceral / Renal		Cerebral		Cerebral		Cerebral		Cerebral		
	J-Curve 1 (RC1)		J-Curve 2 (RC2)		Bentson 1 (JB1)		Mani (MAN)		Newton 3 (HN3)		Newton 4 (HN4)		
Diameter	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	
5 Pack	5 F	100	503110202468	100	503110202469	100	503110202511	100	503110202507	100	503110202503	100	503110202504
	6 F	100	503110202768	100	503110202769	100	503110202811	100	503110202807	100	503110202803	100	503110202804
	7 F	100	503110203068	100	503110203069	100	503110203111	100	503110203107	100	503110203103	100	503110203104
	8 F	100	503110203368	100	503110203369	100	503110203411	100	503110203407	100	503110203403	100	503110203404

	Visceral / Renal		Visceral / Renal		Visceral / Renal		Visceral / Renal		Visceral / Renal		Visceral / Renal		
	Renal 55 (R55)		Renal 65 (R65)		Renal DC Right (RDC A1R)		Renal DC Left (RDC A1L)		Renal DC Right (RDC A2R)		Renal DC Left (RDC A2L)		
Diameter	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	
5 Pack	5 F	100	503110202435	100	503110202434	100	503110202454	100	503110202455	100	503110202456	100	503110202457
		65		65		65	503110202529	65	503110202530	65	503110202531	65	503110202532
	6 F	100	503110202735	100	503110202734	100	503110202754	100	503110202755	100	503110202756	100	503110202757
		65		65		65	503110202829	65	503110202830	65	503110202831	65	503110202832
	7 F	100	503110203035	100	503110203034	100	503110203054	100	503110203055	100	503110203056	100	503110203057
		65		65		65	503110203129	65	503110203130	65	503110203131	65	503110203132
	8 F	100	503110203335	100	503110203334	100	503110203354	100	503110203355	100	503110203356	100	503110203357
		65		65		65	503110203429	65	503110203430	65	503110203431	65	503110203432

	Visceral / Renal		Visceral / Renal		Uterine		Visceral / Renal		Visceral / Renal		Cerebral		
	Shepard Hook 0.8 (SHK 0.8)		Shepard Hook 1.0 (SHK 1.0)		BK		Chueng B (CHG B)		Chueng C (CHG C)		CK 1 (CK 1)		
Diameter	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	Length	Catalog Number	
5 Pack	5 F	100	503110202451	100	503110202450	100	503110202482	100	503110202476	100	503110202477	100	503110202480
		60	503110202520	60	503110202522								
		45	503110202519	45	503110202521								
	6 F	100	503110202751	100	503110202750	100	503110202782	100	503110202776	100	503110202777	100	503110202780
		60	503110202820	60	503110202822								
		45	503110202819	45	503110202821								
	7 F	100	503110203051	100	503110203050	100	503110203082	100	503110203076	100	503110203077	100	503110203080
		60	503110203120	60	503110203122								
		45	503110203119	45	503110203121								
	8 F	100	503110203351	100	503110203350	100	503110203382	100	503110203376	100	503110203377	100	503110203380
		60	503110203420	60	503110203422								
		45	503110203419	45	503110203421								

www.alvimedica.com

Please contact your local sales representative for product availability.

Not for sale in the USA and Japan.

The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.

©2014, Alvimedica. All Rights Reserved.



back to life!