


Alvipro™ 418


0,018" OVER THE WIRE PTA BALLOON CATHETER

The Alvipro™ 418 OTW PTA Balloon Catheter is designed to let you manage the most challenging endovascular procedures.



- Excellent crossability and rewrap allow for multiple inflations under difficult circumstances
- 4 Fr. compatibility for the whole range to minimize the puncture site
- Superior abrasion and scratch resistance for a remarkable performance in calcified vessels

TECHNICAL DATA

BALLOON	
Diameters:	2 - 7 mm
Lengths:	15 - 220 mm
Balloon material:	Polyamide 12
Markers:	2 Platinum/Iridium
Tip material:	Pebax
Tip Entry Profile:	0.026"
DEVICE SPECIFICATIONS	
Over The Wire (OTW)	
Dual Lume	
Compliance	Semi-Compliant
Introducer sheath compatibility:	4F
Maximum guidewire compatibility:	0.018"
Usable Catheter Length - UCL:	90 & 150 (Ø 2-3.5 mm) 90 & 135 (Ø 4-7 mm)
Nominal Pressure:	6 atm
Rated Burst Pressure:	11 - 18 atm

PRESSURE CHART

L (mm) / Ø (mm)	15	20	30	40	60	80	100	120	150	220
2.0		6/18	6/18	6/18	6/18	6/18	6/18	6/18	6/18	6/18
2.5		6/18	6/18	6/18	6/18	6/18	6/18	6/18	6/18	6/18
3.0		6/16	6/16	6/16	6/16	6/16	6/16	6/16	6/16	6/16
3.5		6/16	6/16	6/16	6/16	6/16	6/16	6/16	6/16	6/16
4.0	6/16	6/16	6/16	6/16	6/16	6/16	6/16	6/16	6/16	6/16
5.0	6/14	6/14	6/14	6/14	6/14	6/14	6/14	6/14	6/14	6/14
6.0		6/12		6/12	6/12	6/12	6/12	6/12	6/12	6/12
7.0		6/11		6/11	6/11	6/11	6/11	6/11	6/11	

Nominal Pressure/Rated Burst Pressure

Alvipro™ 418

0,018" OVER THE WIRE PTA BALLOON CATHETER

ORDER INFORMATION ALVIPRO™ 418

USABLE CATHETER LENGTH 90 CM


∅ (mm) \ L	15	20	30	40	60	80	100	120	150	220
2.0		512110206401	512110206402	512110206403	512110206404	512110206405	512110206406	512110206407	512110206408	512110206409
2.5		512110206410	512110206411	512110206412	512110206413	512110206414	512110206415	512110206416	512110206417	512110206418
3.0		512110206419	512110206420	512110206421	512110206422	512110206423	512110206424	512110206425	512110206426	512110206427
3.5		512110206428	512110206429	512110206430	512110206431	512110206432	512110206433	512110206434	512110206435	512110206436
4.0	512110206437	512110206438	512110206439	512110206440	512110206441	512110206442	512110206443	512110206444	512110206445	512110206446
5.0	512110206447	512110206448	512110206449	512110206450	512110206451	512110206452	512110206453	512110206454	512110206455	512110206456
6.0		512110206457		512110206458	512110206459	512110206460	512110206461	512110206462	512110206463	512110206464
7.0		512110206465		512110206466	512110206467	512110206468	512110206469	512110206470	512110206471	

USABLE CATHETER LENGTH 150 CM

∅ (mm) \ L (mm)	15	20	30	40	60	80	100	120	150	220
2.0		512110206501	512110206502	512110206503	512110206504	512110206505	512110206506	512110206507	512110206508	512110206509
2.5		512110206510	512110206511	512110206512	512110206513	512110206514	512110206515	512110206516	512110206517	512110206518
3.0		512110206519	512110206520	512110206521	512110206522	512110206523	512110206524	512110206525	512110206526	512110206527
3.5		512110206528	512110206529	512110206530	512110206531	512110206532	512110206533	512110206534	512110206535	512110206536

USABLE CATHETER LENGTH 135 CM

∅ (mm) \ L (mm)	15	20	30	40	60	80	100	120	150	220
4.0	512110206601	512110206602	512110206603	512110206604	512110206605	512110206606	512110206607	512110206608	512110206609	512110206610
5.0	512110206611	512110206612	512110206613	512110206614	512110206615	512110206616	512110206617	512110206618	512110206619	512110206620
6.0		512110206621		512110206622	512110206623	512110206624	512110206625	512110206626	512110206627	512110206628
7.0		512110206629		512110206630	512110206631	512110206632	512110206633	512110206634	512110206635	



- 1 - Central lumen (guide wire lumen)
- 2 - Balloon lumen
- 3 - Balloon (nominal length)
- 4 - Marker rings
- 5 - Tip
- 6 - Shaft, working length

www.alvimedica.com

Please contact your local sales representative for product availability.

Not for sale in the USA and Japan.

The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.

©2014, Alvimedica. All Rights Reserved.



back to life!