


Alvipro™ 535


0,035" OVER THE WIRE PTA BALLOON CATHETER

The Alvipro™ 535 OTW PTA Balloon Catheter is a semi-compliant workhorse balloon that enables you to be successful in the majority of everyday cases.



- Excellent crossability and rewrap allow for multiple inflations under difficult circumstances
- A-traumatic tapered tip for smooth trackability
- Superior abrasion and scratch resistance for a remarkable performance in calcified vessels

TECHNICAL DATA

BALLOON	
Diameters:	3 - 10 mm
Lengths:	20 - 150 mm
Balloon material:	Polyamide 12
Markers:	2 Platinum/Iridium
Tip Material:	Pebax
Tip Entry Profile:	0.040"
DEVICE SPECIFICATIONS	
Over The Wire (OTW)	
Dual Lume	
Compliance	Semi-Compliant
Introducer sheath compatibility:	5F & 6F
Maximum guidewire compatibility:	0.035"
Usable Catheter Length - UCL:	75 & 120 cm
Nominal Pressure:	6 atm
Rated Burst Pressure:	9 - 16 atm

PRESSURE CHART

∅ (mm) \ L (mm)	20	40	60	80	100	120	150
3,0	6/16	6/16	6/16	6/16	6/16	6/16	6/16
4,0	6/16	6/16	6/16	6/16	6/16	6/16	6/16
5,0	6/14	6/14	6/14	6/14	6/14	6/14	6/14
6,0	6/12	6/12	6/12	6/12	6/12	6/12	6/12
7,0	6/12	6/12	6/12	6/12	6/12	6/12	6/12
8,0	6/11	6/11	6/11	6/11	6/11*	6/11*	6/11*
9,0	6/9*	6/9*	6/9*	6/9*	6/9*	6/9*	6/9*
10,0	6/9*	6/9*	6/9*	6/9*	6/9*		

Nominal Pressure/Rated Burst Pressure * 6Fr Compatible

Alvipro™ 535

0,035" OVER THE WIRE PTA BALLOON CATHETER


ORDER INFORMATION ALVIPRO™ 535

USABLE CATHETER LENGTH 75 CM

∅ (mm) \ L (mm)	20	40	60	80	100	120	150
3,0	512110206801	512110206802	512110206803	512110206804	512110206805	512110206806	512110206807
4,0	512110206808	512110206809	512110206810	512110206811	512110206812	512110206813	512110206814
5,0	512110206815	512110206816	512110206817	512110206818	512110206819	512110206820	512110206821
6,0	512110206822	512110206823	512110206824	512110206825	512110206826	512110206827	512110206828
7,0	512110206829	512110206830	512110206831	512110206832	512110206833	512110206834	512110206835
8,0	512110206836	512110206837	512110206838	512110206839	512110206840	512110206841	512110206842
9,0	512110206843	512110206844	512110206845	512110206846	512110206847	512110206848	512110206849
10,0	512110206850	512110206851	512110206852	512110206853	512110206854		

USABLE CATHETER LENGTH 120 CM

∅ (mm)	20	40	60	80	100	120	150
3,0	512110206901	512110206902	512110206903	512110206904	512110206905	512110206906	512110206907
4,0	512110206908	512110206909	512110206910	512110206911	512110206912	512110206913	512110206914
5,0	512110206915	512110206916	512110206917	512110206918	512110206919	512110206920	512110206921
6,0	512110206922	512110206923	512110206924	512110206925	512110206926	512110206927	512110206928
7,0	512110206929	512110206930	512110206931	512110206932	512110206933	512110206934	512110206935
8,0	512110206936	512110206937	512110206938	512110206939	512110206940	512110206941	512110206942
9,0	512110206943	512110206944	512110206945	512110206946	512110206947	512110206948	512110206949
10,0	512110206950	512110206951	512110206952	512110206953	512110206954		



- 1 - Central lumen (guide wire lumen)
- 2 - Balloon lumen
- 3 - Balloon (nominal length)
- 4 - Marker rings
- 5 - Tip
- 6 - Shaft, working length

www.alvimedica.com

Please contact your local sales representative for product availability.

Not for sale in the USA and Japan.

The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.

©2014, Alvimedica. All Rights Reserved.



back to life!