



Isthmus Logic

CoCr BALLOON EXPANDABLE STENT FOR ILIAC VESSELS



Isthmus Logic is a Cobalt-Chromium stent with Bio Inducer Surface technology for the treatment of iliac vessels.



- Proven bio and haemo-compatibility enabled by the Bio Inducer Surface coating, which accelerates the rate of endothelization and reduces thrombogenicity
- Uniform radial force and excellent scaffolding with an unique cell architecture
- Platinum markers and no stent foreshortening allow unmatched placement accuracy and visualization

TECHNICAL DATA

PERIPHERAL STENT	
Stent diameters:	6 - 10 mm
Stent lengths:	19 - 59 mm
Stent material:	L605 CoCr alloy
Strut surface coating:	Bio Inducer Surface
2 platinum radiopaque markers at stent ends	
No stent foreshortening upon expansion	
DEVICE SPECIFICATIONS	
Over The Wire (OTW)	
Introducer sheath compatibility:	6F
Guidewire compatibility:	0.035"
Usable Catheter Length - UCL:	75 & 135 cm
Nominal Pressure:	7 - 8 atm
Rated Burst Pressure:	10 - 12 atm

Isthmus Logic

CoCr BALLOON EXPANDABLE STENT FOR ILIAC VESSELS


ORDER INFORMATION ISTHMUS LOGIC

USABLE CATHETER LENGTH 75 CM

∅ (mm) \ L (mm)	19	29	39	49	59
6.0	ICLC6019S	ICLC6029S	ICLC6039S		
7.0	ICLC7019S	ICLC7029S	ICLC7039S	ICLC7049S	ICLC7059S
8.0	ICLC8019S	ICLC8029S	ICLC8039S	ICLC8049S	ICLC8059S
9.0	ICLC9019S	ICLC9029S	ICLC9039S	ICLC9049S	ICLC9059S
10.0	ICLC10019S	ICLC10029S	ICLC10039S		

USABLE CATHETER LENGTH 135 CM

∅ (mm) \ L (mm)	19	29	39	49	59
6.0	ICLC6019L	ICLC6029L	ICLC6039L		
7.0	ICLC7019L	ICLC7029L	ICLC7039L	ICLC7049L	ICLC7059L
8.0	ICLC8019L	ICLC8029L	ICLC8039L	ICLC8049L	ICLC8059L
9.0	ICLC9019L	ICLC9029L	ICLC9039L	ICLC9049L	ICLC9059L
10.0	ICLC10019L	ICLC10029L	ICLC10039L		



www.alvimedica.com

Please contact your local sales representative for product availability.

Not for sale in the USA and Japan.

The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.

Manufactured by: CID S.p.a., Strada per Crescentino, sn - 13040 Saluggia (VC) Italy

©2014, Alvimedica. All Rights Reserved.



back to life!